

IHI RECOMMENDATIONS

- “For the operator placing the central line and for those assisting in the procedure, maximal barrier precautions means strict compliance with hand hygiene and wearing a cap, mask, sterile gown, and gloves. The cap should cover all hair and the mask should cover the nose and mouth tightly. These precautions are the same as for any other surgical procedure that carries a risk of infection.
- For the patient, applying maximal barrier precautions means covering the patient from head to toe with a sterile drape, with a small opening for the site of insertion.”

- Institute for Healthcare Improvement,
100,000 Lives Campaign,
“Prevent Central Line Infections: How-to Guide”

CDC GUIDELINES

- “Maximal sterile barrier precautions (e.g., cap, mask, sterile gown, sterile gloves, and large sterile drape) during the insertion of CVCs substantially reduces the incidence of CRBSI compared with standard precautions...”

- Centers for Disease Control and Prevention
“Guidelines for the Prevention of Intravascular
Catheter-Related Infections,” CDC Recommendations
and Reports, August 9, 2002 / 51(RR10);1-26

JCAHO REQUIREMENTS

- “Medications and solutions both on and off the sterile field should be labeled even if there is only one medication being used.”

- Joint Commission on Accreditation of Healthcare Organizations,
2006 Patient Safety Goals

Vantex is a trademark of Edwards Lifesciences Corporation.

BIOPATCH is a registered trademark of ETHICON, INC.

SorbaView, Centurion and SnagFree are registered trademarks of Centurion Medical Products Corporation.

Oligon is a trademark of Implemed, Inc.

ChloraPrep is a registered trademark of Enturia®.

CENTURION
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843
517.546.5400 | 800.248.4058 | fax 517.546.9388 | www.centurionmp.com

©2009 Centurion Medical Products Corporation LIT020V1

CENTURION®




CVC INFECTION REDUCTION BUNDLE

Designed to provide the practitioner everything needed to complete central line insertions

The Centurion® CVC Infection Reduction Bundle doesn't just give you the components necessary to minimize infection risk, it also lets you choose those components that best meet your needs.

Featuring:

- Edwards Vantex Antimicrobial Catheters
- Centurion SorbaView® Dressings
- The BIOPATCH® Protective Disk with CHG
- Centurion Catheter Securement System
- Centurion Catheter Anchor Block
- Centurion SnagFree® Needle Holder
- ChloraPrep®
- Drapes, gowns, masks and other components necessary for optimum protection against bloodstream related infections



VANTEX CATHETER WITH OLIGON

The Vantex catheter features Oligon which provides a Broad Spectrum Antimicrobial Protection.

The Vantex has an infusion rate up to 77% faster than competing catheters.



CENTURION CVC SECUREMENT ANCHOR

Centurion CVC Securement Anchor allows for easy application.




CENTURION CVC BUNDLE

The Centurion CVC Bundle organizes all of your components in an easy to use system. Drapes, gowns and caps are easily configured over the tray compartment.

HOSPITAL FORMS

Your hospital forms can be built into your custom trays.



SALINE AND HEPARIN FLUSHES

A large selection is available.



* Pre-filled syringes are attached to outer layer.

CHLORAPREP 3ML APPLICATOR





SAFETY COMPONENTS

A large selection of safety scalpels and needles from major manufacturers are available.



VALVES

Safety valves and closed valve systems from major manufacturers are available.



SNAGFREE® NEEDLE HOLDER

The Centurion SnagFree Needle Holder has a patented design to prevent the suture from snagging on the instrument while tying knots. SnagFree Instruments are put through a rigorous quality control process and are guaranteed to equal the performance of more expensive, reusable substitutes.

Also available with large loops.

SORBAVIEW® ULTIMATE IJ DRESSING SYSTEM

The SorbaView Ultimate IJ Dressing System with anchor is the one product that you can be certain will not fall off your patient. Its specially formulated adhesive provides a proven barrier to liquids and bacteria, while allowing excellent breathability. SorbaView is easy to apply, minimizes patient discomfort, and remains intact for extended periods.




SORBAVIEW ULTIMATE WITH BIOPATCH DISK


Our newest version of the SorbaView Ultimate Dressing has an adhesive-free zone that is designed to fit over the BIOPATCH Protective Disk. This system increases the area of protection for the insertion site and allows for easy removal of the dressing without jeopardizing the integrity of the BIOPATCH.



SORBAVIEW 2000

The SorbaView 2000 is designed specifically for insertions at subclavian and chest sites. Its specially formulated adhesive provides a proven barrier to liquids and bacteria, while allowing excellent breathability and minimizing patient discomfort.





DRAPES FOR CVC INSERTIONS

In addition to offering a large variety of stock drapes, Centurion can fully customize a drape to suit your specific needs.

A Size: 40" x 65"
 Material: Nonwoven
 Fenestration: 4" with adhesive border
 Perforation: Yes
 Absorbent Zone: 32" x 33"


B Size: 65" x 104"
 Material: Nonwoven
 Fenestration: 3/4" with adhesive border
 Perforation: Yes
 Window: 18" x 18" clear plastic

C Size: 60" x 90"
 Material: SMS
 Fenestration: 3" with adhesive border
 Window: 10" x 10" clear plastic

D Size: 65" x 90"
 Material: Nonwoven
 Fenestration: 4" with adhesive zone


E Size: 60" x 90"
 Material: SMS
 Fenestration: 4"
 Window: 8" circle adhesive zone

Material Key:



PROTECTING CAREGIVERS

Centurion can supply your every need, from catheters to gloves, you can rely on the quality of Centurion components and products.



GOWNS

- Surgical gowns
- Multiple protection levels available
- Non-reinforced, fabric reinforced and impervious
- Multiple sizes available



STERILE LABELS

Sterile labels can be custom made to fit your exact procedural requirements.

