

Protect light cord with convenient pre-attached sheath in 3 quick steps!



Insert light source into handle



Loosen rubber band from sheath



Pull sheath down over cord



CENTURION®

# PELVIC EXAM TRAYS

AND DISPOSABLE

# VAGINAL SPECULA



**CENTURION**  
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843  
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | [www.centurionmp.com](http://www.centurionmp.com)

0626-09

**CENTURION**  
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843  
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | [www.centurionmp.com](http://www.centurionmp.com)

