

CENTURION®


SORBAVIEW® 2000 WINDOW DRESSING

DESIGNED WITH NURSING RECOMMENDATIONS TO...

GO ON EASIER

LAST LONGER

REDUCE COST OUTCOMES



(includes tape strip and a HubGuard®)


SV30NVST (includes tape strip)

See your Centurion® sales representative for product demonstrations and trials.

A comprehensive line of Centurion catheter securement products are available to further help stabilize and protect sites.



INCLUDE SORBAVIEW® 2000
IN YOUR CUSTOM TRAY



¹Lee, M. Seven Day Dwell Times for Vascular Access Device Exit Site Dressings, Lee Medical Inc. 2005.

²Penney-Timmons, E. Decreased Costs/Improved Outcomes With Standardized Intravenous Equipment. *Journal of the Association for Vascular Access (JAVA)*, Vol 10, No 1 (2005), pp 20-23.

³Garcia, R et al. NY hospital uses process, product education to cut bloodstream infections. *Healthcare Purchasing News.*, Sept. 2003.

⁴Garcia, R et al. Three Years Experience in Implementing HICPAC Recommendations for the Reduction of Central Venous Catheter-related Bloodstream Infections. Brookdale University Medical Center (BUMC), 2003.

U.S. Patent No. 6,841,715

CENTURION
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | www.centurionmp.com

LIT-SV2-500

0610-08


CENTURION®

SORBAVIEW® 2000 WINDOW DRESSING

CENTURION®

SORBAVIEW® 2000 WINDOW DRESSING

- Easy to apply – won't stick to itself during application
- Stays securely in place
- Provides proven barrier against bacteria and contaminants
- Absorbent pad wicks away moisture
- One-way permeable film allows skin to breathe
- Increases patient comfort²
- Comprehensive offering of sizes
- Non-toxic, non-sensitizing, non-irritating and latex free
- Available in individual sterile packages, or in custom procedure trays designed to help prevent phlebitis, infiltration and reduce risk of infection^{2,3,4}



SorbaView® 2000 Window Dressing sets the standard of excellence in catheter site management and securement.


SorbaView® is easy to apply, and remains intact for up to seven days, reducing dressing changes and thus lowering the risk of infection.^{1,3,4}

SorbaView® allows skin to breathe, while providing a proven barrier to liquids and bacteria. Its signature design allows vapors to escape through the film window while the absorbent pad wicks away excess moisture. This patented action prevents build-up of fluid at the site – a breeding ground for bacteria – and keeps moisture from breaking down the adhesive under the fabric border. Adhesive stays dry, dressings stay on!

Available in a large assortment of sizes for managing a variety of catheter sites.

WHAT MAKES SORBAVIEW® DRESSINGS BETTER...

OUR CONSTRUCTION



- Meets CDC recommendations for a seven day dressing change protocol¹
- Centurion® Dressings protect the site placing the occlusive film layer directly on the skin's surface – contaminants cannot penetrate the barrier!
- One-way permeable film exceeds normal skin respiration rates, ensuring that skin breathes normally under the dressing.
- Other dressings place porous fabric next to the skin, potentially compromising the dressing's seal and creating a pathway to the site!