


THE CENTURION® PICC SECURITY SYSTEM™



THREE POINTS OF SECURITY
IN ONE COMPLETE SYSTEM

1. stabilizes the catheter
2. protects the patient
3. secures the site

CENTURION
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | www.centurionmp.com

THE CENTURION[®] PICC SECURITY SYSTEM

A complete, all inclusive system for PICC Securement

1. stabilizes the catheter

- Adhesive-based securement device stabilizes virtually any size PICC
- Pliable construction conforms comfortably to body contours
- Low profile design improves dressing adherence

2. protects the patient

- Provides proven bacterial barrier
- Remains intact for up to seven days, reducing dressing changes and thus lowering the risk of infection

3. secures the site

- Dressing specifically sized to fit the securement device
- One-way permeable film allows skin to breathe while the absorbent pad wicks away excess moisture
- Tape strip helps seal and secure the site more effectively



Security System Advantages

- Simple application and easy removal – dressing won't stick to itself during application or to securement device upon removal
- Standardizes procedures – takes the guesswork out of application and removal
- Non-sensitizing, non-irritating, latex free
- Meets INS Standards for Catheter Stabilization
- U.S. Patent No. 6,841,715, additional patent pending