



CENTURION®

# Suction Specimen Traps

## Baby-friendly Catheter

- Catheter tips are formed with round edges and hand-polished to be soft, smooth and non-abrasive
- Two small recessed eyelets near the tip improve suctioning and reduce clogging
- Catheters have smooth satin finish to reduce insertion trauma

## Clinician-friendly Hydrophobic Filter

- Filter prevents contamination of suction canisters and tubing, reducing supply costs and time spent changing canisters
- Trap functions on its side, even if mucus contacts part of the filter
- Reduces the risk of exposure to body fluids during emergency oral suctioning

## Specimen Trap

- Container seals tightly for safe transport to lab
- Trap is graduated to 20cc

*Safely obtain, measure and transport  
meconium or mucus specimens*

**CENTURION**  
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843  
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | [www.centurionmp.com](http://www.centurionmp.com)

# CENTURION®



## Mucus Trap with Catheter

- For use with wall suction, portable suction, or emergency oral suctioning
- Basic sterile set includes pre-assembled device and label, gloves are optional
- Available with 6, 8, 10, 12 or 14 Fr. catheter



## Meconium Aspirator

- Attaches to any endotracheal tube
- Thumb port for regulating suction
- Packaged sterile



## Dual Function Meconium Aspirator / Mucus Trap

- Instead of stocking two products in the Delivery Room, now only one is needed!
- As packaged, this sterile device is used for normal suctioning and specimen collection
- If intubation is required, simply separate at the adapter and attach meconium aspirator to endotracheal tube for suctioning

**CENTURION**  
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843  
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | [www.centurionmp.com](http://www.centurionmp.com)

0625-09